Mini-Plus

User's Manual

Care for a Healthy Life

CONTENTS

COMPRESSOR NEBULIZER

USER'S MANUAL 1

COMPRESOR NEBULIZADOR

MANUAL DE INSTRUCCIONES 17

COMPRESSOR NEBULIZADOR

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MODEL NO.: 9R-021 SERIES

Please read the manual before use

CAUTION - U.S. Federal Law restricts this device to sale by or on the order of a physician.

This product has been tested and successfully approved to the following standards:

EN 60601-1

EN 60601-1-2

0197

IMPORTANT SAFEGUARDS SAVE THESE INSTRUCTIONS

READ ALL INSTRUCTIONS BEFORE USING

DANGER - To reduce the risk of electrocution:

1. Always unplug this product immediately after using.
2. Do not use while bathing, showering, dish washing, or close to water sources of any kind.
3. Do not place or store product where it can fall or be pulled into a tub or sink.
4. Do not place in or drop into water or other liquid.
5. Do not reach for a product that has fallen into water. Unplug immediately.

WARNING - To reduce the risk of burns, electrocution, fire or injury to persons:

1. This product should never be left unattended when plugged in.
2. Close supervision is necessary when this product is used by, on or near children or invalids. Choking accident may result from a child swallowing a small part that has become detached from the device or its accessories.
3. Use this product only for its intended use as described in this manual. Use this product only under doctor's direction. Do not use attachments not recommended by the manufacturer.
4. Never operate this product if a) it has a damaged cord or plug, b) it is not working properly, c) it has been dropped or damaged, or d) it has been dropped into water. Return the product to a specified service center for examination and repair.

5. Do not wrap the power cord around the compressor (main unit).
6. Disconnect the power plug by pulling the plug, not by pulling on the compressor (main unit), or the cord.
7. If the power cord or plug becomes frayed or otherwise damaged, do not use.
8. Do not place heavy objects on the power cord, or bend and pull the cord harder than necessary. These actions could cause an electric shock or fire.
9. Keep the cord away from heated surfaces.
10. Never block the air openings of this product or allow objects to fall or be inserted into the air vent openings or place it on a soft surface such as bed or couch, where the air openings may be blocked.
11. Never use while sleeping or feeling drowsy.
12. Never drop or insert any object into any opening or hose.
13. No modification of this equipment is allowed.
14. Do not use in outdoors or operate where aerosol (spray) products are being used or where oxygen is being administered in a closed environment such as an oxygen reservoir.
15. Potential allergic reactions to accessible materials used in the Compressor Nebulizer equipment. If any signs of allergic reaction or hypersensitivity happen, stop the treatment immediately, and notify the doctor or nurse.
16. Potential contact injuries for patients used in the Compressor Nebulizer equipment. If any contact injuries happen, stop the treatment, and notify the doctor or nurse.
17. Do not leave long lengths of air tubing around the top of your bed. It could twist around your head or neck while you are sleeping and lead to strangulation.
18. When using any electrical appliance certain important safety measures must always be observed, including the following:

NOTE, CAUTION, WARNING STATEMENTS AND SYMBOLS:

NOTE - Indicates information that user should pay special attention to.

CAUTION - Indicates correct operating or maintenance procedures in order to prevent damage to or destruction of the equipment or other property.

WARNING - Calls attention to a potential danger that requires correct procedures or practices in order to prevent personal injury.

EMC Warning Statement

This equipment has been tested and found to comply with the limits for medical devices to the IEC 60601-1-2. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving device.
- Increase the separation between the equipment.
- Connect the equipment into an outlet on a circuit different from that to which the other device(s) are connected.
- Consult the manufacturer or field service technician for help.

SYMBOLS

	Authorised representative in the European community
	Manufacturer
	Off, disconnection from the mains
	On, connection to the mains
IP21	Protected against solid foreign objects of 12,5 mm and greater; Protection against vertically falling water drops (Only for 230V System)
	Alternating current
	Temperature Limitation
	CAUTION should be used to highlight the fact that there are specific warnings or precautions associated with the device
	"BF" symbol, indicate this product is according to the degree of protecting against electric shock for type BF equipment
	Class II
	Consult operating instructions for use.
	Consult instructions for use
	Disposal of Electrical & Electronic Equipment (WEEE): This product should be handed over to an applicable collection point for the recycling of electrical and electronic equipment. For more detailed information about the recycling of this product, please contact your local city office, household waste disposal service or the retail store where you purchased this product.

1. Introduction

1.1 Intended Use

The Mini Plus Compressor Nebulizer System is intended to provide a source of compressed air for aerosol therapy. It is used in conjunction with a jet (pneumatic) nebulizer to produce medicated aerosols for inhalation by pediatric and adult patients with respiratory symptoms.

CAUTION - Indications for therapy include asthma, chronic bronchitis, infection of the upper respiratory tract, chronic obstructive pulmonary disease and other respiratory disorders in accordance with a medical doctor's prescription. Other than the usage mentioned above, please do not use this product for any other purpose. This device can be used with adults or pediatric patients under physician's prescription.

1.2 Safety Precaution Instruction

When using this electrical product, especially when children are present, one should always follow basic safety precautions. Do not install, maintain or operate this equipment without reading, understanding and following the proper Mini Plus Compressor System instruction manual, otherwise injury or damage may result. Important information is highlighted by using the following bulletin:


DANGER : Urgent safety information for hazards that will cause serious injury or death.

WARNING : Important safety information for hazards that might cause serious injury or property damage.

CAUTION : Information for preventing damage to the product.

2. Product Description

1. Power Switch
2. Filter Cap (disposable filter inside)
3. Air-Outlet Connector
4. Integrated Carrying Handle
5. Power Cord
6. Tubing Holder
7. Cooling Air Openings

3. Operation

NOTE - Before first use, the nebulizer cup assembly should be cleaned following instructions described in the "Cleaning" section.

WARNING - Before connecting the power cord, make sure the I/O (ON/OFF) switch is in the 0 (OFF) position.

The plug is also used to disconnect the device. Do not position the equipment so that it is difficult to disconnect the device.


The device is also equipped with a plug to enable connection to the electric power supply.

3.1 Daily Use Operation

CAUTION - The Mini Plus Compressor is designed for intermittent use only. Do not operate it continuously for more than 30 minutes for a single use without turning it off and allowing a cooling period of at least 30 minutes.

1. Before each use inspect the Mini Plus Compressor and nebulizer cup assembly for damage or wear, and replace as needed.
2. Place the Mini Plus Compressor on a table or other flat stable surface. Be sure

you can easily reach the controls when seated. Do not use this device on the floor.

3. With the power switch in the 0 (OFF) position, plug the power cord into an appropriate electrical wall outlet.
4. Connect one end of the tubing to the compressor air-outlet connector.
5. Assemble the nebulizer cup and add indicated medication to the nebulizer cup before use.
6. Attach the other end of the tubing into the air-inlet connector at the bottom of the nebulizer cup.
7. Turn on the Compressor by pressing the power switch to the I (ON) position and begin treatment.
8. If treatment needs to be interrupted, simply press the power switch to the 0 (OFF) position.
9. When the treatment is complete, turn off the compressor by pressing the power switch to the 0 (OFF) position and unplug the unit from the electrical outlet.

Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide. This device has no AP or APG protection.

3.2 Safety Overload

NOTE - Do not exceed 30 minutes of continuous operation.

1. The motor of this device has a built-in thermal overload protector. Should the motor overheat, the protector will automatically shut off the motor. Should this occur, turn the I/O (ON/OFF) switch to the 0 (OFF) position and allow the motor to cool down for approximately 30 minutes before turning it on again.

2. If the overload protector shuts off the motor frequently, you may have an unstable voltage situation.
3. If the unit shuts down and cannot restart, it may need to be replaced. Call your provider immediately.

4. Cleaning

4.1 Compressor Outer Case Cleaning

WARNING - Electric shock hazard. Do not remove outer case of this unit. All disassembly and maintenance of this unit must be done by a qualified service technician. Refer servicing to qualified service personnel.

WARNING - This unit does not require oil. Do not attempt to lubricate any internal parts.

WARNING - Unplug unit before cleaning. Do not submerge in water for cleaning.

1. Wipe the main unit with a damp cloth every few days to keep it dust-free.
2. Do not use any powdered cleaners or soap. Do not submerge the unit into water.

4.2 Nebulizer Cup Cleaning

All parts of the disposable nebulizer, except tubing, should be thoroughly cleaned after each use. If your doctor specifies a specific cleaning procedure, follow that procedure.

WARNING - To prevent possible risk of infection from contaminated medication, cleaning of the nebulizer is recommended after each treatment.

1. After each use, disassemble the nebulizer cup, tube and mouthpiece completely, wash all parts in warm soapy water and rinse thoroughly. Allow to air dry. Store in a clean, dry location.

2. Once daily, after washing as above, soak nebulizer components in a fresh solution consisting of one part warm water and one part white vinegar for 30 minutes followed by a warm water rinse.

NOTE - The disposable nebulizer kit is for single patient use only.

NOTE - Please follow national requirements to dispose of the unit properly.

5. Storage

Keep the unit and accessories dry. Avoid direct sunlight.

6. Maintenance



6.1 General Information

1. This unit is oil-less. Do Not Lubricate.
2. Risk of electric shock. Do not disassemble the main unit.
3. All maintenance must be performed by qualified service personnel.

6.2 Filter Change

NOTE - Change the filter once every month or when filter turns gray. Please follow the instructions below as shown to the right.

1. Open the filter cap.
2. Remove dirty filter with a small, pointed object. Discard the filter.
3. Replace with a clean filter. Additional filters should be

purchased from your provider.

4. Put the filter cap back.

CAUTION - Do not use cotton or any other material as a filter. Do not wash or reuse the filter. Do not operate the unit without a filter.

7. Troubleshooting

If your Mini Plus Compressor Nebulizer fails to function, consult the Troubleshooting guide below. If the problem persists, consult your equipment provider.

Problem	Cause and Solution
Device doesn't operate.	<ul style="list-style-type: none"> • Check if plug is properly fitted into an appropriate electrical outlet. • When device has been running continuously for over 30 minutes right before using, an auto shut down may be activated by the built-in thermal protector; cool down device for 30 minutes before next usage.
Weak Nebulization	<ul style="list-style-type: none"> • Check for proper electrical voltage. • Check tubing for blockage or air leakage at connection to Mini Plus compressor or nebulizer cup, and replace as needed. • Check the nebulizer cup to see if it is properly assembled and not damaged. If there is any damage, replace as needed. • Check if filter is too dirty, replace as needed.

8. Specifications

Electrical Rating (Note: Refers to the rating label on the product)	<ul style="list-style-type: none"> ◆ 120VAC,60Hz,1.2A ◆ 230VAC,60Hz,0.6A 	<ul style="list-style-type: none"> ◆ 230VAC,50Hz,0.6A
Maximum Compressor Pressure	≥ 35 psi	≥ 35 psi
Operating Flow Rate	≥ 6.0 lpm	≥ 5.5 lpm
Sound Level	≤ 72 dBA at 1 meter distance (Device only)	≤ 70 dBA at 1 meter distance (Device only)
Aerosol Output	≥ 1.5 ml	
Aerosol Output Rate	≥ 0.3 ml/min	
Medicine cup capacity	Minimum 2ml - Maximum 6ml	
Particle size	0.5 ~ 5 microns*	
Classification	Class II. BF equipment. IP21. No AP/APG protection.	
Applied part	Mouthpiece or Mask	
Dimensions (W x D x H)	14.0 × 18.8 × 10.0 cm / 5.5" × 7.4" × 3.9"	
Weight	1.6 kg / 3.5 lb	
Fuse (see note 3)	F5AH 250V for 120V system T1.6AL 250VAC for 230V system	
Expected Service Life	3 Year	
Environment	Temperature	Operation: 10°C to 40°C / 50°F to 104°F
		Storage: -15°C to 50°C / 5°F to 122°F
		Transport: -15°C to 70°C / 5°F to 158°F
	Humidity	Operation: 10% to 90%RH non-condensing
		Storage: 10% to 90%RH non-condensing
		Transport: 10% to 90% RH non-condensing
Atmospheric pressure	Operation: 700-1013.25 hPa	

- Note:**
1. *Values are determined by supplier's nebulizer.
 2. Consult distributor(s) or EU representative(s) for additional technical documents.
 3. The Nebulizer is fitted with a fuse as specified in the Specifications chart, which is designed to be maintained by the manufacturer or its agents. Users are not allowed to open the cover to replace a blown fuse.

9. Accessories

<u>Model</u>	<u>Discription</u>
668064-0000	Jet Nebulizer Set (Adult & Child Masks)
668024-0000	Jet Nebulizer Set (Mouthpiece)

APPENDIX A: EMC INFORMATION

Guidance and Manufacturer's Declaration- Electromagnetic Emissions:

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment-Guidance
RF emissions CISPR 11	Group1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network
Harmonic emissions IEC61000-3-2	Class A	
Voltage fluctuations / Flicker emissions IEC61000-3-3	Complies	

Warning:

1. The device should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device should be observed to verify normal operation in the configuration in which it will be used.
2. Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or

decreased electromagnetic immunity of this equipment and result in improper operation.
 3. Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Nebulizer, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Guidance and Manufacturer’s Declaration- Electromagnetic Immunity:

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Basic EMC standard	Immunity Test Levels		Compliance Levels	Electromagnetic Environment-Guidance
	Professional healthcare facility environment	HOME HEALTHCARE ENVIRONMENT		
Electrostatic Discharge (ESD) IEC61000-4-2	±8kV contact ±15kV air		±8kV contact ±15kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/ burst IEC61000-4-4	±2kV for power supply line ±1kV for input/output line		±2kV for power supply line ±1kV for input/output line	Mains power quality should be that of atypical commercial or hospital environment
Surge IEC61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s)	± 1 kV line(s) to line(s)	Mains power quality should be that of atypical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	Voltage Dips: i) 100% reduction for 0.5 period, ii) 100% reduction for 1 period, iii) 30% reduction for 25/30 period, Voltage Interruptions: 100% reduction for 250/300 period		120V / 230V	Mains power quality should be that of a typical commercial or hospital environment. If the user of this device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC61000-4-8	30 A/m	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

<p>Conducted RF IEC 61000-4-6</p>	<p>3 Vrms 0,15 MHz – 80 MHz 6 Vrms in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz</p>	<p>3 Vrms 0,15 MHz – 80 MHz 6 Vrms in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz</p>	<p>6Vrms</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of this device, including cables, than there commended separation distance calculated from the equation applicable to the frequency of the transmitter.</p>
<p>Radiated RF EM Fields IEC61000-4-3</p>	<p>3 V/m 80 MHz to 2.7 GHz 80 % AM at 1 kHz 385-6000 MHz, 9-28V/m, 80% AM(1kHz) pulse mode and other modulation</p>	<p>10 V/m 80 MHz to 2,7 GHz 80 % AM at 1 kHz 385-6000 MHz, 9-28V/m, 80% AM(1kHz) pulse mode and other modulation</p>	<p>10V/m</p>	<p>Recommended separation distance $d = \sqrt{P}$ 150kHz to 80MHz $d = 0.6\sqrt{P}$ 80MHz to 800MHz $d = 1.2\sqrt{P}$ 800 MHz to 2.7G MHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).^b Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey , should be less than the compliance level in each frequency ranged. Interference may occur in the vicinity of equipment marked with the following symbol: </p>

NOTE 1: U_T is the a.c. mains voltage prior to the application of the test level

NOTE 2: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 3: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people

- a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.
- b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

Recommended separation distances between portable and mobile RF communications equipment and this device:

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = \sqrt{P}$	80 MHz to 800 MHz $d = 0.6\sqrt{P}$	800 MHz to 2,7 GHz $d = 1.2\sqrt{P}$
0.01	0.1	0.06	0.12
0.1	0.31	0.19	0.38
1	1	0.6	1.2
10	3.1	1.9	3.8
100	10	6	12

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

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Print-2017/All rights reserved