# XT Fit

## CPAP System

### *Instruction Manual*

#### USA

APEX MEDICAL USA CORP.  
615 N. Berry St. Suite D, Brea, CA 92821, USA

APEX MEDICAL S.L.  
Elcano 9, 6ª planta  
48008 Bilbao. Vizcaya. SPAIN

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## IMPORTANT SAFEGUARDS – PLEASE MAKE THIS PRECAUTIONS LIST AVAILABLE TO ALL USERS and PROVIDERS

READ ALL INSTRUCTIONS BEFORE USING THIS SYSTEM

### WARNINGS:

- THIS DEVICE IS NOT INTENDED FOR LIFE SUPPORT.
- If oxygen is used with the device, the oxygen flow must be turned off when the device is not in use.  
**Explanation of the Warning:** When the device is not in operation and the oxygen flow is left on, oxygen delivered into the tubing may accumulate within the device's enclosure and create a risk of fire.
- Oxygen supports combustion. Oxygen should not be used while smoking or in the presence of an open flame.
- Always ensure the device is generating airflow before the oxygen supply is turned. Always turn off the oxygen supply before stopping the airflow from the device.
- This device **should not** be used in the vicinity of a flammable anesthetic mixture in combination with oxygen or air and nitrous oxide.
- The airflow for breathing generated by this device may be as much as 7°C (12.6° F) higher than the room temperature. This device **should not** be used if the room temperature is warmer than 35°C (95°F) to prevent the airflow temperature from exceeding 40° C (104° F) and causing irritation to your airway.
- If this device overheats, it will stop operating and show “**Error 002**” on the display. After cooling down to proper temperature, the device can restart again.
- This machine should be used only with masks (and connectors) recommended by the manufacturer, or by your physician or respiratory therapist. A mask should not be used unless the CPAP machine is turned on and operating properly. The vent holes associated with the mask should never be blocked for proper exhaling purpose. If the vent hole is blocked, the CPAP machine will stop and show message “**Error 002**”, after cooling down, please re-connect the power cord to reset the machine.
- At low CPAP pressure, some exhaled gas may remain in the mask and be re-breathed.

### CAUTION –

- Make sure the environment around the machine is dry and clean. Dust and foreign particles may affect the treatment. Keep the air inlet on the back of the machine clear to prevent overheating and damage of the device. Do not place the machine near a source of hot or cold air. Extreme cold or hot environment may damage user's respiratory airway.
- If there is a possibility of electro-magnetic interference with mobile phones, please increase the distance between devices or turn off the mobile phone.
- Do not connect the device to the personal computer for data downloading during the treatment. This may cause the CPAP system failure.
- To prevent from potential allergic reactions to mask used in the device, This machine shall be used only with masks (and connectors) recommended by the manufacturer.
- U.S. Federal law restricts this device to sale by or on the order of a licensed physician.

**DANGER** –Reducing the risk of electrocution:

1. Always unplug this product immediately after using.
2. **DO NOT** use while bathing.
3. **DO NOT** place or store product where it can fall into or be pulled into a tub or sink.
4. **DO NOT** place or immerse the unit into water or other liquids.
5. **DO NOT** reach for the unit if it is immersed in water. Unplug the unit immediately.

**WARNINGS** - To reduce the risk of burns, electrocution, fire or injury:

1. This product should never be left unattended when plugged in.
2. Close supervision is necessary when this product is used by, on, or near children or invalids. Choking resulting from a child swallowing a small part that has become detached from the device or its accessories.
3. Use this product only for its intended purpose as described in this manual. **DO NOT** use auxiliary attachments that are not approved by the manufacturer.
4. Never operate this product if it has a damaged cord or plug. If it is not working properly, or has been dropped, damaged, or immersed in water, return the product to a service center for examination and repair.
5. Keep the cord away from hot surfaces and objects.
6. Never block the air outlets on this product. **DO NOT** place the unit on a soft surface such as a bed or couch, where the air outlets may become obstructed. Keep the air outlets free of lint, hair, and other obstructive particles.
7. Never drop or insert any object into any opening on the unit or into the tubing.
8. When disposing the unit because of damage or failure, you must comply with the regulations of your particular region or country for disposing such equipment.
9. For therapy pressure 4-18 cmH<sub>2</sub>O, the default pressure compensation is available for altitude from sea level up to 8000 feet; for therapy pressure 18.5-20 cmH<sub>2</sub>O, the default pressure compensation is available for altitude from sea level up to 5000 feet.
10. No modification of this equipment is allowed.
11. Do not modify this equipment without authorization of the manufacturer.
12. If this equipment is modified, appropriate inspection and testing must be conducted to ensure continued safe use of the equipment.
13. Do not leave long lengths of air tubing around the top of your bed. It could twist around your head or neck while you are sleeping and lead to strangulation.

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## 1. Introduction

This manual is provided to explain initial set up and operating instructions. This manual should be stored in a safe location and be readily accessible for future reference.

### 1.1 General Information

Obstructive Sleep Apnea (OSA) is a condition that is an intermittent and reoccurring obstruction of the upper respiratory tract, causing a complete (apnea) or partial (hypopnea) blockage of respiratory airflow during sleep. The syndrome varies in degree, depending on the position of the tongue and soft palate muscle when in a relaxed state.

The most common treatment for OSA is Continuous Positive Airway Pressure (CPAP). CPAP devices can deliver constant air pressure into your upper respiratory tract via a nasal mask. This constant air pressure can keep your airway open during sleep, preventing occurrences of OSA.

This unit is a CPAP device regulated by an internal micro-processor. It features a menu-driven LCD display, universal power supply, and Ramp Time Adjustment. The Ramp Time Adjustment feature allows the unit to begin at low pressures, allowing you to comfortably fall asleep. The unit will then gradually "ramp up" the pressure over a configurable period of time until the appropriate therapy pressure is reached. This unit also has a Compliance Meter that records patient usage data that can be reviewed by a physician at a later time.

The system has been tested and successfully approved for the following standard:

0197

EN 60601-1  
EN 60601-1-2  
EN 61000-3-2 Class A  
EN 61000-3-3

**For US and CANADA only**

E348970  
53DG

Medical Equipment- CPAP  
**WITH RESPECT TO ELECTRICAL SHOCK,  
FIRE AND MECHANICAL HAZARDS ONLY**  
IN ACCORDANCE WITH ANSI/AAMI  
ES60601-1 (2005,3rd ed.) and CAN / CSA C22.2  
No. 60601-1 (2008)

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## EMC Warning Statement

This equipment has been tested and found to comply with the limits for medical devices for **EN 60601-1-2:2007**. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. This unit generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the unit.
- Increase the distance between the conflicting equipment.
- Connect the unit into an outlet that is on a different circuit than the other device(s).
- Consult the manufacturer or field service technician for help.

This system has been tested to comply with the following voluntary standards: **FDA**

## 1.2 Intended Use

This device is intended to provide continuous positive airway pressure (CPAP) for the treatment of adult Obstructive Sleep Apnea (OSA).

**Cautions:** Some patients might have pre-existing contraindications for CPAP therapy, or might experience some potential side effects of using CPAP device, please consult your physician if you have any questions concerning your therapy.

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## 2. Product Description

Components including:

1. Main CPAP device
2. Detachable power cord
3. User manual
4. Flexible air tubing with 1.8 m length
5. Carrying bag

### 3. Installation

#### 3.1 Unpacking

To secure its contents inside, the CPAP device and accessories are bundled in a paper-packaged box. Unpack this box by removing the CPAP and its accessory and checking for any damage, which may have occurred during shipping. If there are damages, please contact your dealer immediately.

#### 3.2 Setting Up

- 1) Connect the power cord to CPAP device and plug into main electrical outlet.

Once the power cord is plugged into the electrical outlet, the device is in ready to operate position ("STANDBY" sign appears in LCD display)

**NOTE:** The plug is also served to disconnect the device.

- 2) Connect one end of the air tubing firmly onto the air outlet of the CPAP.

- 3) Connect the other end of the air tubing to the mask system. Putting on the mask and headgear according to the mask instruction manual.

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## 4. Operation

**⚠️ NOTE:** Always read the operating instruction before use.

### 4.1 Control Panel Description

Main control panel and button descriptions:

#### START/STANDBY

To start the treatment, simply press the **"START/STANDBY"** button. To stop the treatment, press the **"START/STANDBY"** button again. The display will switch between **[STANDBY]** and Therapy Pressure **[XX.X cmH<sub>2</sub>O]** in cmH<sub>2</sub>O unit.

#### MENU

When the unit is in "Standby" mode, you can press the **"MENU"** button to enter the configuration menus. The patient can adjust the "Ramp Up Time" (Ramp), Ramp Starting Pressure (Ramp P) and Altitude Compensation (Alt). The three other options, Therapy Pressure (P), Compliance Meter (CM), and Total Operating meter (TM) are to be set only by a physician or other care provider. To confirm a change to a setting, you must press the **"MENU"** button again. Refer to section 4.2 (Menu Option Descriptions) for a detailed explanation.

#### UP

Press the **"UP"** button to the increasing values when adjusting configurations.

#### DOWN

Press the **"DOWN"** button to decrease values when adjusting configurations.

## 4.2 Getting Started – First time use.

For first-time use, follow the steps below:

1. Make sure that your physician or care provider has initialized and configured the appropriate treatment values into the unit.
2. Make sure that your package contains a face or nasal mask with headgear. If not, your physician or care provider can recommend the type of mask and headgear you must obtain
3. Attach tubing and mask (see Setup section).
4. Attach power cord to the unit and an electrical outlet (see Setup section). Unit will power-on automatically.
5. Use the **"MENU"** button to set the Current Time, Ramp Up Time, Ramp Starting Pressure, and Altitude Compensation (see **"MENU"** Button Section for description of these options). Note that these values may have been already set by your physician or care provider.
6. Press the **"START/STANDBY"** button – Airflow to the mask will begin.
7. You can now place the Headgear and Mask over your head and face and begin treatment.

## 4.3 Menu Button - Option Descriptions

The following items describe the options presented to you after you press the **"MENU"** button on the control panel.

### (1) Ramp Time (Ramp XX MIN)

The Ramp Time feature on this unit allows the patient to fall sleep at a lower, more comfortable pressure, after which the unit will gradually ramp up to the appropriate therapy pressure. The first option of the **"MENU"** button is the Ramp Time, which is the time (in minutes) it will take for the unit to "ramp" up to the appropriate therapy pressure. Press the **"UP"** or **"DOWN"** buttons to set the preferred ramp time and press **"MENU"** for confirmation. There are ten timing options; 0, 5, 10, 15, 20, 25, 30, 35, 40 and 45 minutes.

### (2) Ramp Starting Pressure (Ramp P XX.X)

Press **"MENU"** button to select [Ramp P XX.X] menu, press **"UP"** or **"DOWN"** button to set the preferred ramp starting pressure and press **"MENU"** for confirmation. The ramp starting pressure can be changed from 3 cmH<sub>2</sub>O to "Therapy Pressure – 1" cmH<sub>2</sub>O. For example, if your therapy pressure is 10 cmH<sub>2</sub>O, the maximum ramp starting pressure you can select is 9 cmH<sub>2</sub>O.

### (3) Therapy Pressure (P XX.XcmH<sub>2</sub>O)

The third option of the **"MENU"** button is the Therapy pressure displayed in units of cmH<sub>2</sub>O. The Therapy pressure option should only be set by the equipment provider, a respiratory therapist or a physician.

 **NOTE: The Therapy Pressure value itself should only be prescribed by a physician.**

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#### (4) Altitude Compensation (Alt XXX)

The fourth option of the **"MENU"** button is the Altitude Compensation option, used to adjust pressure in the mask based on the altitude relative to sea level. Press the **"UP"** or **"DOWN"** buttons to select from eight options (1 to 8) that correspond to a range of altitudes listed in the table below. The level should be set relative to your elevation above sea level. Once the preferred option has been selected, press the **"MENU"** button to confirm.

NOTE: Users can operate XT Fit at wide range of altitude within 0~8000 ft (limit to 5000 ft if the pressure is set at 18.5-20 cmH<sub>2</sub>O). The altitude function provides a method of accuracy pressure output at high altitude. Users must follow the below table to select correct altitude setting depending on your elevation above sea level, the device will automatically regulate airflow output to achieve targeted pressure at high altitude. Otherwise, wrong altitude setting will cause inaccuracy pressure output.

There are eight altitude compensation levels. The level should be set depending on your elevation above sea level. Below table is the comparison of altitude setting and altitude.

Altitude Setting	Altitude (Imperial)	Altitude (Metric)
1	0 ~ 1000 ft	0 ~ 304 m
2	1001 ~ 2000 ft	305 ~ 609 m
3	2001 ~ 3000 ft	610 m ~914m
4	3001 ~ 4000 ft	915 m ~1219m
5	4001 ~ 5000 ft	1220 m ~1524m
6	5001 ~ 6000 ft	1525 m ~1829m
7	6001 ~ 7000 ft	1830 m ~2134m
8	7001 ~ 8000 ft	2135 m ~2438m

#### (5) Total Meter (TM XXXX.X)

The fifth **"MENU"** item is the Total Meter which displays the total number of hours that the device has been "On" and delivering airflow. In other words, it is the time the **"START/STANDBY"** button was pressed to begin airflow to the time the **"START/STANDBY"** button was pressed again to stop airflow. This meter should only be re-set by the equipment provider, a respiratory therapist or by a physician.

#### (6) Compliance Meter (CM XXXX.X)

Press **"MENU"** button to select **[CM XXXX.X hr]** menu, the compliance meter records the total therapy hours for the device. The compliance meter should be re-set only by the provider, a respiratory therapist or by a physician.

#### (7) Session Counter (CM XXX 4hr)

Press **"MENU"** button to select **[CM XXX 4hr]** menu, the session counter displays the number of sessions greater than four hours. Your home care provider may ask for this information. The session counter will be erased once the compliance meter has been reset.

**(8) Check Value (c XXX)**

Press **"MENU"** button to select [**c XXX**] menu, which displays a check value. Your home care provider may ask for this information.


**(9) Turn off the Device**

Remove the power cord from the electrical outlet, and disconnect power cord from the power socket on the back of device.

**⚠ NOTE:** Once the setting is confirmed, press **"MENU"** button. Otherwise, the device will automatically go back to standby without saving the modification if no action is taken in 20 seconds.

**4.4 Flowchart of Menu settings**

Enter the user's menu mode by pressing the **"MENU"** button.

In each setting, when the preferred value has been selected, press **"MENU"** for confirmation and press **"MENU"** again to enter next selection.

**⚠ NOTE:** For physicians, please refer to a separated "Physician's Additional Instruction" page.

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## 5. Adding a Humidifier

XT Fit CPAP device can be used with XT Heated Humidifier **(9S-006500)** which is available from the home care provider. The heated humidifier may reduce nasal dryness and irritation by providing adequate moisture and heat to the airflow. Please refer to the XT heated humidifier **(9S-006500)** instruction manual for complete setup information.

**NOTE:** When XT Fit CPAP device is used with the heated humidifier, its power supply is from the power socket outlet of the heated humidifier. **DO NOT** connect the power cord of CPAP device and plug into main electrical outlet.

## 6. Cleaning & Maintenance

### 6.1 Cleaning the Main Unit

The main unit should be cleaned regularly (at least every 30 days). To clean the unit, wipe it down with a damp cloth using a mild detergent. When selecting a cleaning agent, take precautions to choose one that is not chemically abrasive to the surface of the plastic case. All parts should be air-dried thoroughly before use. During storage, the unit should be kept in a dust-free location.

**WARNING:** Don't try to open this device. Repairs and internal servicing should only be performed by an authorized service agent. Don't drop any subject into the air tube or air outlet.

### 6.2 Tubing and Mask

The tubing and mask should be checked and cleaned daily. Please refer to the cleaning instruction packaged with the accessories.

1. Disconnect the air tubing from the air outlet of the device.
2. Remove the air tubing and headgear straps from the nasal mask.
3. Use a mild detergent and prepare the detergent according to manufacturer's recommendations. Tubing and mask can be washed and rinsed in warm water.
4. Rinse thoroughly and make sure they are completely dry before next use.
5. All items of the mask and air tubing are subject to normal wear and tear and may eventually be replaced. Replace the mask and the air tubing if they are damaged.

**CAUTION:** Do not use bleach, chlorine-, alcohol-, or aromatic-based (including all scented oils), moisturizing or antibacterial soaps to clean the cushion, mask, air tubing. These solutions may cause hardening and reduce the life of the product.

- CAUTION:** Do not wash or dry the mask or air tubing at a temperature above 70°C (160°F)
- WARNING:** Do not use any cleaner containing fragrance or conditioners as they will leave a residue.
- WARNING:** The mask must not be re-used by another person. This is to avoid the risk of cross-infection.

### 6.3 Air Filter Maintenance

For an optimum operation of the device, the air filter can be cleaned by the user at least once every two weeks or more often if this device is operated in a dusty environment. It is recommended for the filter to be replaced with a new one after 6 months of use of the device; however, depending on the air environment, the replacement time may vary.

- CAUTION:** Dirty air filter may cause high operating temperatures that affect device performance. Ensure the air filter is cleaned and fitted at all times.

1. Remove the dirty filter from the enclosure on the rear of the device.
2. Wash the filter in warm water with a mild detergent, and rinse with water. Allow the filter to air dry completely before reinstalling. **DO NOT** use a filter that is not completely dry. If the filter is torn, replace it.
3. Reinstall the filter.

- NOTE:** Please follow national requirements to dispose the unit properly.

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## 7. Troubleshooting

The table below lists troubleshooting solutions for the problems that may happen. If the problem persists, contact your equipment provider service agent.

Problem	Possible Causes	Solutions
No display	<ol style="list-style-type: none"> <li>1. The power cord is not connected to the power socket.</li> <li>2. LCD failure or controlled PCB failure.</li> </ol>	<ol style="list-style-type: none"> <li>1. Ensure the power cord is connected.</li> <li>2. Contact your equipment provider for repair.</li> </ol>
Display code incorrect	LCD failure or controlled PCB failure.	Contact your equipment provider for repair.
Illuminant under LCD is not on	LED failure	Contact your equipment provider for repair.
Buttons disable	Button failure	Contact your equipment provider for repair.
Air delivered is slow	<ol style="list-style-type: none"> <li>1. During ramp time.</li> <li>2. Filter is too dirty.</li> <li>3. Flow generator failure.</li> </ol>	<ol style="list-style-type: none"> <li>1. Check the ramp time setting</li> <li>2. Change or clean the filter regularly.</li> <li>3. Contact your equipment provider for repair.</li> </ol>

### Error / Warning Messages show in LCD

Message type	Definition	Message in LCD
Error: Primary function can't execute.	Error for abnormal system settings	<b>Error 001</b>
	Error for flow generator failure	<b>Error 002</b>
Warning:	Out of system memory	<b>Warn 001</b>
	<ol style="list-style-type: none"> <li>1. System memory is nearly full</li> <li>2. Memory resetting problem</li> </ol>	<b>Warn 002</b>

**⚠** **NOTE:** When the warning message appears, contact your physician or equipment provider to download the memory data and reset the meter.

## 8. Technical Specifications

Item	Specifications	
Mode of Operation	Continuous	
Pressure Range	4 –20 cmH <sub>2</sub> O (adjustable in 0.5 cmH <sub>2</sub> O increment)	
Ramp Time	0 – 45 minutes (adjustable in 5-minute increment)	
Ramp Starting Pressure	3 to therapy pressure	
Altitude Compensation (Manual setting)	Level 1~8 for 0 ~ 8000 ft when the pressure set at 4-18 cmH <sub>2</sub> O. Level 1~5 for 0 ~ 5000 ft when the pressure set at 18.5-20 cmH <sub>2</sub> O	
Dimensions (W x D x H)	14.5 x 13.0 x 10.0 cm or 5.7" x 5.1" x 3.9"	
Weight	Approximately 800 g or 1.76 lb	
Sound Level	28 dBA at 10 cmH <sub>2</sub> O, 1 meter distance	
Power Input	AC100-240V, 50/60 Hz, 0.5-0.3A	
Environment	Temperature	Operating: +5° C to +35° C (+41° F to +95° F) Storage/Transport: -15° C to +50° C (+5° F to +122° F)
	Humidity	Operating: 15%RH to 95%RH non-condensing Storage/Transport: 10%RH to 90%RH non-condensing
Air Tubing	Flexible plastic, 1.8m (approx.)	
Classification:	Class II Type BF, Applied Parts Nasal Mask. Not suitable for use in the presence of a flammable anesthetic mixture. IP21: Protected against solid foreign objects of 12.5 mm and greater; Protected against vertically falling water drops. Continuous operation. Continuous operation.	

**Note:** The manufacturer reserves the right to modify the specification without notice.

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**Symbols**

BF symbol, which indicated this product, is according to the degree of protecting against electric shock for type BF equipment.

Attention, should read the instructions.

Refer to instruction manual

Attention, should read the instructions.

Class II

**IP21** 

Protected against solid foreign objects of 12.5 mm and greater; Protected against vertically falling water drops.

Disposal of Electrical & Electronic Equipment (WEEE):

This product should be handed over to an applicable collection point for the recycling of electrical and electronic equipment. For more detailed information about the recycling of this product, please contact your local city office, household waste disposal service or the retail store where you purchased this product.

Authorized representative in the European community

Manufacturer

**9. NOTE, CAUTION, AND WARNING STATEMENTS**

**NOTE:** Indicate information that you should pay special attention to.

**CAUTION:** Indicate correct operating or maintenance procedures in order to prevent damage to or destruction of the equipment or other property.

**WARNING:** Calls attention to a potential danger that requires correct procedures or practices in order to prevent personal injury.

## 10. Servicing

The XT series systems are intended to provide safe and reliable operation according to the instructions provided by Apex Medical. Apex Medical recommends that the XT system be inspected and serviced by authorized technicians if there is any sign of wear or concern with device function. Otherwise, service and inspection of the devices generally should not be required during the five year design life of the device.

## 11. Limited Warranty

Apex Medical Corp. (hereafter 'Apex') gives the customer a limited manufacturer warranty on new original Apex products and any replacement part fitted by Apex in accordance with the warranty conditions applicable to the product in question and in accordance with the warranty periods from date of purchase as listed below. This warranty does not cover damage caused by accident, misuse, abuse, alteration and other defects not related to material or workmanship. To exercise your rights under this warranty, please contact your local, authorized Apex dealer.

Product	Warranty period
CPAP devices	2 years
Water chamber	6 months
Mask systems	6 months
Disposable products	None

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## Appendix A: EMC Information

### Guidance and Manufacturer's Declaration- Electromagnetic Emissions:

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment-Guidance
Harmonic emissions IEC61000-3-2	Class A	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network.
Voltage fluctuations / Flicker emissions IEC61000-3-3	Complies	

### Guidance and Manufacturer's Declaration- Electromagnetic Immunity:

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Immunity Test	IEC60601 test level Compliance		Electromagnetic Environment-Guidance
Electrostatic Discharge(ESD) IEC61000-4-2	±6kV contact ±8kV air	±6kV contact ±8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC61000-4-4	±2kV for power supply line ±1kV for input/output line	±2kV for power supply line ±1kV for input/output line	Mains power quality should be that of atypical commercial or hospital environment
Surge IEC61000-4-5	±1kV for differential mode ±2kV for common mode	±1kV for differential mode ±2kV for common mode	Mains power quality should be that of atypical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	<5 % UT(>95 % dip in UT)for 0,5 cycle 40 % UT(60 % dip in UT)for 5 cycles 70 % UT(30 % dip in UT)for 25 cycles <5 % UT(>95 % dip in UT)for 5 sec	<5 % UT(>95 % dip in UT) for 0,5 cycle 40 % UT(60 % dip in UT) for 5 cycles 70 % UT(30 % dip in UT) for 25 cycles <5 % UT(>95 % dip in UT) for 5 sec	Mains power quality should be that of atypical commercial or hospital environment. If the user of this device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of atypical location in a typical commercial or hospital environment.

NOTE:  $U_T$  is the a.c. mains voltage prior to the application of the test level

**Guidance and Manufacturer's Declaration- Electromagnetic Immunity:**

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Immunity Test	IEC60601 test level	Compliance	Electromagnetic Environment-Guidance
Conducted RF IEC 61000-4-6	3Vrms150 kHz to 80 MHz outside ISM bands <sup>a</sup>	3Vrms	Portable and mobile RF communications equipment should be used no closer to any part of this device, including cables, than there commended separation distance calculated from the equation applicable to the frequency of the transmitter.  Recommended separation distance $d = 1.2\sqrt{P}$ 150kHz to 80MHz
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3V/m	$d = 1.2\sqrt{P}$ 150kHz to 80MHz  $d = 2.3\sqrt{P}$ 80 MHz to 2.5GHz  Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). <sup>b</sup>  Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>c</sup> , should be less than the compliance level in each frequency range <sup>d</sup> .  Interference may occur in the vicinity of equipment marked with the following symbol:  

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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- a) The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz.
- b) The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2.5 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 is used in calculating the recommended separation distance for transmitters in these frequency ranges.
- c) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.
- d) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

**Recommended separation distances between portable and mobile RF communications equipment and this device**

This device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this device as recommended below, according to the maximum output power of the communications equipment

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1,2\sqrt{P}$	80 MHz to 800 MHz $d = 1,2\sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance *d* in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**Note 1:** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**Note 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

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Print-2012/All rights reserved